

COMPUTERISED 1 & 3 MHZ ULTRASOUND THERAPY



Advanced and trend setter Ultrasound therapy from HMS R&D. DIGISONIC 3s is designed indigenously with international standards with incomparable quality. The DIGISONIC 3S comes with 1 MHZ and 3 MHZ outputs with contact control sensor facility. Micro controller based circuitry with LCD display for viewing all parameters. 47 preprograms and 30 re-writable user settable memory. Contact controls suspends the Ultrasound output when transducer contact with skin becomes insufficient. Easily selectable continuous and 9 pulsed modes soft touch keys for all functions.

FEATURES

1 MHZ and 3 MHZ outputs in single transducer

Contact control sensor

Graphic LCD display

Evidence based Pre-programmes

User programmable memory

Water resistant transducer

Energy Efficient

Two size transducer

Compact and stylish

Feather touch Cap-sense keys

EN ISO 13485: 2012 + AC: 2012

SPECIFICATIONS

Frequency - 1MHZand3MHZ
Therapymodes - Continuous/Pulsed

Pulse mode - 10,20,30,40,50,60,70,80,

90&100%

PulseRate - 100Hz

Intensity - $0.1-3 \, \text{W/cm}^2$

ERATransducer L - 5cm²

ERATransducerS - 0.8cm²

Pre-programmes - 47 Evidence based protocals User memory - 30 Re-writable positions

LCD Display - 128x64 pixels

Timer - 1-99 mins

TECHNICAL DATA

Mains Voltage - 90-270V/50 or 60 Hz

Safety Class - Class I, Type BF (IEC 60601 -1)
Dimensions - L x B x H - 23 x 22.8 x 9.5 cms

Weight - 1.7 Kgs

Web: www.hms.co.in email: sales@hms.co.in