

DIGITENS 104



DIGITENS 104 is an advanced micro controller based Four Channel TENS from HMS R & D. Unique design with manual and 20 pre-programs. Used for both fitness and therapeutic purpose. Soft touch keys for selection of programs. Built in battery with 10 hrs backup indigenous design with international standards.

FEATURES


- 4 Channel Outputs
- Micro controller based circuitry
- Manual Cum preprograms facility
- Obesity Program
- Universal Input
- Built In Re-chargeable battery with 10 hrs Back-up
- ABS Cabinet-Portable & Light Weight

TECHNICAL DATA

- Mains Input Voltage - 90V – 270V AC
- Battery - Maintenance free rechargeable
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H – 21 x 22.5 x 6.5 cms
- Weight - 1.75 Kgs

SPECIFICATIONS

- Tens Frequency - 2 – 250Hz (1 Hz Step)
- Pulse width - 20 – 250 micro seconds
- No.of Channels - Four
- Output mode - Continuous/Pulsed
- Frequency Sweep - Available
- Treatment Time - 0 – 60 min
- Pre-program - 1-15 Pain Relieving Program
16-20 Fitness Program
- Display - 128 / 64 Graphic LCD Display
- Power Source - Dual Power (Mains/Battery)



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in