

## HYDROCOLLATOR



### FEATURES

Durable Stainless steel water tank  
Thermostatically Controlled  
Three size of silica gel packs  
Mobile unit fitted with wheel  
Simple to fill and drain water

### SPECIFICATIONS

Tank Capacity - 15 Litres  
Heater type - Industrial heater  
Tank Made up of - Stainless steel  
Watts - 1500 w

### TECHNICAL DATA

Mains Voltage - 220V/50Hz  
Dimensions - LxBxH - 36x23x26 cms  
Weight - 19 Kgs (without water)