

DIGISTIM EMS

MICRO STIMULATOR - POCKET MODEL



DIGISTIM EMS is a specially designed Pocket model Stimulator unit for physiotherapy students and home visit therapists. It has three modes IG, F & SF. IG duration / SF intervals can be varied by single control. Micro-controller based circuitry helps us to make it compact and sleek. It can be operated by battery/mains.

FEATURES

Battery/Mains operated Microcontroller based circuitry Handy model

SPECIFICATIONS

Treatment Models - IG, F, SF
IG - 30–150 mSec.

FARADIC

Freq./P.width - 50 Hz / 0.7 mSec

Surge Duration - 1-6 sec

Power Source - 9 V Battery / Adaptor



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116. Ph: +91 44 2476 0756, 2476 6576 Telefax: +91 44 2476 3773

Web: www.hms.co.in email: sales@hms.co.in