

**COMPUTERISED 1 & 3 MHZ ULTRASOUND THERAPY**



Advanced and trend setter Ultra sound therapy from HMS R&D. DIGISONIC is designed indigenously with international standards with incomparable quality. The DIGISONIC comes with 1 MHZ and 3 MHZ outputs. Micro controller based circuitry with large LCD display for viewing all parameters. 30 pre-programmes and 20 re-writable user settable memory. Easily selectable continuous and 9 pulsed modes soft touch keys for all functions.

**FEATURES**

- 1 MHZ and 3 MHZ outputs in single transducer
- Large LCD display
- 30 Pre-programmes
- 20 User settable memory
- Water resistant transducer
- User friendly
- Light weight plastic moulded cabinet
- Cap-sense keys
- Plastic moulded water resistant transducer

**SPECIFICATIONS**

- Frequency - 1 MHZ and 3 MHZ
- Therapymodes - Continuous/ Pulsed
- Pulse settings/Hz - 9 Settings/ 100Hz
- Output power - 0-3 W/cm<sup>2</sup>
- Timer - 0-60 mins

**TECHNICAL DATA**

- Mains Voltage - 90-270V/ 50 or 60 Hz
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 28 x 22 x 9 cms
- Weight - 2.3 Kgs

**HMS Medical Systems**

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.  
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773  
Web : [www.hms.co.in](http://www.hms.co.in) email : [sales@hms.co.in](mailto:sales@hms.co.in)