

EXERCISE BANDS



FEATURES

Available with different resistance
Latex and Latex free models available

PHYSIO BALLS



FEATURES

Sizes from 65cm to 120 cm Diameter

GEL BALLS



FEATURES

Three different hardness
Soft, Medium and Hard

ELECTRIC TILT TABLE



FEATURES

Elevated frames to allow proper lifter positioning
Heavy duty motor lift up to 500 Kg
Powder coated frames that resist rusting and scratching
Remote control for up-down adjustments
Angle indicator
Three extra wide movable patient positioning straps with locking swivel casters
Removable footrest



ISO 13485: 2003 + AC: 2009

HMS MEDICAL SYSTEMS

#9, First Cross Street, New Colony,
Porur, Chennai - 600116.

Ph : 91-44-2476 0756, 2476 6576, 2474 0735

Mobile : 09445040858

Telefax : 91-44-2476 3773

web : www.hms.co.in , physiotherapyequipments.net ,

email : sales@hms.co.in , admin@hms.co.in

EXERCISE THERAPY PRODUCTS



ANKLE EXERCISER



FEATURES

Spring load adjustable
Sturdy & durable

SUPINATOR PRONATOR



FEATURES

Supination and pronation in single unit
Made up of MS with powder coated
360 degree angle rotation
Fixing height and diameter of rotation can be adjustable
Friction can be easily adjustable

EQUIPMENT TROLLEY



FEATURES

MS Powder coated
Fitted with caster wheel with break
Detachable

SHOULDER WHEEL



FEATURES

MS powder coated
Friction Adjustable
Height can be adjusted by a single knob.

TRACTION TABLE



FEATURES

Two section traction table with non friction rolling lumbar section
Traction fittings height can be adjustable
MS with powder coated
Cushion top with resin cover

QUADRICEPS CHAIR



FEATURE

MS powder coated
Chrome plated Lever set
With Backrest or Pench type
Ankle adjustable torque arm hold up to 25 Kg weight.