

TABLE TOP DUAL CHANNEL TENS



INDOTENS is an advanced table top dual channel TENS. INDOTENS is a user friendly equipment with three output modes continuous, burst and polarity alternation, independent controls for each channel.

FEATURES


- Independent channels
- Micro-controller based timer user friendly
- Universal input voltage
- 99 min timer
- Separate intensity control for both channels.

TECHNICAL DATA

- Mains Input Voltage - 90V - 270VAC
- Battery - Maintenance free -rechargeable
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 32 x 22 x 9.5 cms
- Weight - 1.75 Kgs

SPECIFICATIONS

- Frequency - 2 - 250Hz (11 Steps)
- Pulse width - 20 - 250 mic sec (cont. Adj)
- No.of Channels - Two
- Output mode - Continuous/Pulsed/Pol.Alt.
- Pulse Freq - 1-10 Hz
- Treatment Time - 0 - 99 min



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in