

## SHORTWAVE DIATHERMY - 500 WATTS



INDOTHERM 500 is a heavy duty shortwave diathermy equipment can be used in hospitals and big clinics. It is mounted in castors to move from place to place. It has plastic moulded multi-position arms set for setting the treatment place without disturbing the patient. Digital timer for setting the treatment time.

### APPLICATIONS

Rapid reduction of wounds  
Rapid reabsorption of haematomas and oedemas  
Powerful stimulations of the peripheral circulation  
Post – traumatic and post operative disorders such  
arthropathy, Contusion, distortion, hametoma  
Other disorders such as ankylopoletic spondylosis,  
bursltis, Coccycodinia, myalgia (acute) and chronic  
sinusitis

### SPECIFICATIONS

Generator Frequency - 27.12 mHz  
Output Power - 500 W (Adj in 5 steps)  
Timer - 0-99 min Digital Timer

### TECHNICAL DATA

Mains Voltage - 220V/ 50Hz  
Dimensions - L x B x H - 35 x 48 x 75 cms  
Weight - 38 Kgs