

## ULTRASOUND THERAPY



**INDOSONIC 103 is a fast and easy start-up device with an extremely fast way to set-up treatment parameters and it gives quick therapeutic results. A compact and light-weight device for easy transportation. The treatment head is made to fit your hand comfortably. A light ring clearly indicates the contact / non contact and effective transfer of energy. Ultrasound therapy is a proven effective treatment method for numerous pathologies.**

### FEATURES

- Provision for combination therapy
- Underwater treatment transducer
- Compact and light weight
- Extremely quick start-up
- Contact control
- 2 Treatment heads possible
- Energy efficient
- Ultra-flexible cable

### SPECIFICATIONS

- Ultrasound Frequency - **1 MHz**
- Ultrasound Output Modes - **Continuous / Pulsed**
- Pulse Duration / Frequency - **1,2,4 & 5 m secs / 100 Hz**
- Intensity - **0 - 3 W / cm<sup>2</sup> (max)**
- Timer - **0-99 mins with set +, set -, start / pause facility**
- Display - **Power output W/cm<sup>2</sup>**

### TECHNICAL DATA

- Mains Voltage - **90 - 270V / 50 or 60 Hz**
- Safety Class - **Class I, Type BF (IEC 60601 -1)**
- Dimensions - **LxBxH - 20x24x8 cms**
- Weight - **1.75 Kgs**

### HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.  
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773  
Web : [www.hms.co.in](http://www.hms.co.in) email : [sales@hms.co.in](mailto:sales@hms.co.in)