


DIGITAL ULTRASOUND THERAPY 3 MHZ



INDOSONIC - D3 is another advanced and trend setter Ultrasound therapy from HMS R&D. INDOSONIC - D3 is designed indigenously with international standards with incomparable quality. The INDOSONIC - D3 comes with 3 MHz output. Fully Micro controller based circuitry with intelligent skin contact control facility. Soft touch keys for all functions.

FEATURES

- 3 MHz output
- Digital display for timer and watts
- Intelligent sensor for skin contact control
- User friendly
- Light weight plastic moulded cabinet
- Ultra feather touch Cap-sense keys
- Plastic moulded water resistant transducer
- Indicator for contact and non contact with skin

SPECIFICATIONS

- Frequency - 1 MHz (or) 3 MHz
- Therapy modes - Continuous / Pulsed
- Pulse settings - 4 Settings (1:9, 2:8, 2:3, 1:1)
- Output power - 0 - 3 W/cm² (15W in continuous 21W in pulsed)
- Timer - 0 - 60 mins

TECHNICAL DATA

- Mains Voltage - 90 - 270V / 50 or 60 Hz
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 23 x 20 x 8 cms
- Weight - 1.8 Kgs

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in