

SHORTWAVE DIATHERMY - 300 WATTS



INDOTHERM 300 W is a portable shortwave diathermy unit for small clinics. It is maintenance free and user friendly equipment. It has a 99 minute digital timer and automatically cuts off after the treatment time.

SPECIFICATIONS

Frequency	- 27.12 MHz
Output Power	- 300 W (Adj in 5 steps)
Timer	- 0-99 mins at the end of treatment time out put cut off automatically and buzzer sound up

ACCESSORIES

1. Mains Cable	- 1 No.
2. Output Cable	- 1 Set
3. Rubber Pad	- 2 Nos.
4. Instructions Manual	- 1 No.

TECHNICAL DATA

Dimensions	- LxBxH – 38.5 x 40 x 30 cms
Weight	- 22 Kgs