


SIX PROGRAMMED DUAL CHANNEL TENS



DIGITENS is an advanced hand held dual Channel Pocket model TENS in the field of Electrotherapy. The microcontroller based circuitry has six preprogrammed output modes that can be selectable by a soft touch key. It has 15 mins built in timer and it will be operated by battery/ mains.

FEATURES

- Six pre programs
- Micro-controller controlled circuitry
- Battery / Adaptor
- Built-in 15 min timer
- Separate intensity control for both channels.



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.

Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773

Web : www.hms.co.in email : sales@hms.co.in