

ELECTRONIC TRACTION



INDOTRAC is an advanced technology solid state circuitry Static / Intermittent traction unit. Moulded Aluminium gear box and metal gears gives long life to the equipment. It has silence smooth and powerful motor to pull the appropriate weight. The treatment time is controlled by digital timer and has indicators for Rest / Hold time. Separate tension control for cervical and pelvic.

FEATURES

Imported sturdy & noise less motor
Moulded aluminium gear box
Patient safety remote
Easy to operate

TECHNICAL DATA

Mains Voltage - 220V/ 50Hz
Safety Class - Class I, Type BF (IEC 60601 -1)
Dimensions - L x B x H - 28 x 31 x 28 cms
Weight - 11 Kgs

SPECIFICATIONS

Treatment Types - Continuous and Intermittent
Traction force - 3 - 60 Kgs (120Kgs with Doubler)
Hold Time - 20, 40, 60, 80 sec
Rest time - 1, 5, 10, 20 sec
Timer - 0 - 99 min. At the end of treatment time Acoustic signal and automatic reduction Of traction force

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in