

## LONGWAVE DIATHERMY



SPOTCURE Longwave Diathermy is computerised user friendly product from HMS R&D. Longwave Diathermy treatment based on the capacitor field method. It can work in heavy voltage fluctuation. Longwave can use as a alternative for Shortwave Diathermy Applications. It is portable and light weight.

### FEATURES

- Microcontroler based
- Large LCD display
- Touch free capsence keys
- Plastic moulded light weight cabinet
- Wide input voltage
- Electrode size 15, 23 & 30mm

### SPECIFICATIONS

- Output Frequency - 1MHz
- Output Power - 1 Watts/100 Ohm. Adj in 15 steps
- Output Mode - Continuous / Pulsed
- Pulse Freq. - 100 Hz (4Steps)
- Display - 128 x 64 Graphic LCD Display
- Timer - 0-60 min

### APPLICATIONS

- Tendinitis
- Low back pain
- Myalgia
- Arthrosis
- Bursitis
- Ligament injuries
- Ankle distortion

### TECHNICAL DATA

- Mains Supply - 90-270V, 50/60 Hz
- Safety Class - Class I, Type BF (IEC 60601-1)
- Dimensions - LxBxH- 30x25x10.5 cms
- Weight - 1.7 Kgs

