

## MUSCLE STIMULATOR



INDOSTIM - DT is a micro controller based Diagnostic / therapeutic muscle Stimulator unit. Visual / Audio indicator for IG mode. IG duration and SF interval varied by separate control. Light weight and compact in design.

### FEATURES

- Continuous galvanic current
- Constant current output
- Interrupted galvanic with variable pulse duration
- Surge variation in seconds
- Light Weight and Compact
- Buzzer for IG pulse

### TECHNICAL DATA


- Mains Voltage - 90-270V /50 or 60 Hz
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 17 x 28.5 x 9 cms
- Weight - 1.6 Kgs

### SPECIFICATIONS

- Current Forms - G / IG / F / SF
- IG Pulse Duration - 0.1, 0.3, .1, .3, 1, 3, 10, 30, 100, 300, 1000 mSec

### FARADIC

- Freq./P.width - 50 Hz / .7 mSec
- Surge Duration - 1-6 sec (Continuously adjustable)
- Output Current - 0-40 mA



### HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.  
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773  
Web : [www.hms.co.in](http://www.hms.co.in) email : [sales@hms.co.in](mailto:sales@hms.co.in)