

MUSCLE STIMULATOR



INDOSTIM - 201 is a micro controller based therapeutic muscle Stimulator unit. Visual / Audio indicator for IG mode. Soft touch key is used to select the treatment modes. IG duration and SF interval varied by same control. Light weight and compact in design.

FEATURES


- Continuous galvanic current
- Constant current output
- Interrupted galvanic with variable pulse duration
- Surge variation in seconds
- Light Weight and Compact
- Buzzer for IG pulse

SPECIFICATIONS

- Current Forms - G / IG / F / SF
- IG Pulse Duration - 30-300 mSec
- FARADIC
- Freq./P.width - 50 Hz / .7 mSec
- Surge Duration - 1-6 sec (Continuously adjustable)
- Output Current - 0-40 mA

TECHNICAL DATA

- Mains Voltage - 90-270V/50 or 60Hz
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 17 x 18 x 7 cms
- Weight - 1.2 Kgs



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in