



## UN WEIGHING SYSTEM - GAIT TRAINER

HMS Rehab GAIT TRAINING SYSTEM is used to reduce pressure of whole body to lower limbs of the patient, while through hanging belts so as to ensure safety. Mainly used for those patients such as arthritis lower limb's disability, pain and spasm caused by nervous system dysfunction, and help patients recover walking function.

### FEATURES

- Dual motor functions
- Remote control functions
- Hanging scale for weight display
- Heavy duty lockable caster wheel
- Treadmill can be attached
- Lift patient up to 200kgs.

### SPECIFICATIONS

Frame	-	Heavy duty metal frame
Height of the main frame	-	225cm
Weighing capacity	-	250kgs
Hanging lifting range	-	0 - 70cm
Lifting mechanism	-	Electric DC operated imported motor
Recommended treadmill	-	Slow medical treadmill
Operating Voltage	-	220V / 50Hz
Dimensions	-	L x B x H - 100 x 110 x 225cm