

DIGITAL POCKET TENS



DIGITENS II is an advanced micro controller based Two Channel TENS from HMS R & D. Used for both fitness and therapeutic purpose. Soft touch keys for selection of programs. Built in battery. Indigenous design with international standards.

FEATURES

- 2 Channel Outputs
- Micro controller based circuitry
- LCD display
- Soft touch keys
- Independent Intensity control
- Universal Input

SPECIFICATIONS

- Tens Frequency - 2 – 150Hz
- Pulse width - 50 – 300 micro seconds
- No.of Channels - Two
- Output mode - Burst, Normal, MRW, SD, BI Pulse,
- Treatment Time - 1 – 30 min
- Display - Graphic LCD Display
- Power Source - Dual Power (Mains/Battery)

TECHNICAL DATA

- Mains Input Voltage - 90V – 270V AC
- Battery - 9V
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H – 21 x 22.5 x 6.5 cms
- Weight - 1.75 Kgs

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in