



HMS HIGH-LOW TILTING TABLE Fully equipped for Rehabilitation and mobility purposes Electric height adjustment, easy to move with castors

FEATURES

- Elevated frames to allow proper lifter positioning
- Heavy duty motor lift up to 250 Kgs
- Powder coated frames that rusting and scratching
- Remote control for up-down adjustments
- Angle indicator
- Three extra wide movable patient positioning straps with locking swivel casters
- Removable footrest

SPECIFICATIONS

- Safe working load with lifting and tilting capacity of 250kgs
- Electric height and tilting operation with hand switch control
- Adjustable angle dual footboards - positive and negative (+15 to -30)
- Worktable supplied with standard model - optional with Economy model
- Supplied with patient handgrips and set of harnesses
- Height range 47cm to 101cm (18-39") Tilt angle from horizontal to 90 degrees vertical
- Lowers to wheelchair height for ease of patient transfer
- Tilt angle inclinometer fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist