

## SHORTWAVE DIATHERMY - P500



INDOTHERM P500 is a heavy duty shortwave diathermy equipment can be used in hospitals and big clinics. It is mounted in castors to move from place to place. It has plastic moulded multi-position arms set for setting the treatment place without disturbing the patient. Digital timer for setting the treatment time.

### APPLICATIONS

Rapid reduction of wounds  
 Rapid reabsorption of haematomas and oedemas  
 Powerful stimulations of the peripheral circulation  
 Post – traumatic and post operative disorders such as arthropathy, Contusion, distortion, hametoma  
 Other disorders such as ankylopoletic spondylosis, bursltis, Coccycodinia, myalgia (acute) and chronic sinusitis

### SPECIFICATIONS

Generator Frequency - 27.12 MHz  
 Pulse Frequency - 50 - 500 Hz (10 steps)  
 Pulse Width - 400  $\mu$  sec  
 Output Power  
     Continuous - 500 W (Adj in 5 steps)  
     Pulse - 1000 W (Adj in 5 steps)  
 Timer - 0-99 min Digital Timer

### TECHNICAL DATA

Mains Voltage - 220V/ 50Hz  
 Dimensions - L x B x H - 35 x 48 x 75 cms  
 Weight - 38 Kgs