

COMPUTERISED LASER THERAPY



ACCESSORIES



The DIGILASER 203n is an advanced and modern designed unit for laser therapy with interchangeable laser probe with auto detection. Fully computerised 5" Colour LCD with touch screen display. Laser therapy, also known as photo biomodulation, Low Level Laser Therapy (LLLT), cold laser therapy, and as laser biostimulation is an effective method for:

- Ø Pain Management
- Ø Soft Tissue Injury Treatment
- Ø Wound Healing
- Ø Skin Treatment (esp. for Acne)
- Ø Nerve Regeneration
- Ø Cosmetic Regeneration
- Ø Resolving Viral and Bacterial Infections
- Ø Diabetic wound healing

FEATURES

Autoprobe detection
 IR with visible Red laser combination
 Microcontroller based
 5" Colour LCD with touch screen
 Capacitive Touch Display
 Lightweight plastic cabinet

APPLICATIONS

Dermatological Conditions
 Rheumatic Conditions
 Degenerative Joint Disease
 Post-Traumatic and Post-Operative Conditions
 Rheumatoid Arthritis
 Conditions for which acupuncture is indicated

SPECIFICATIONS

Laser type - Semiconductor laser/GaAlAs
 Laser probe
 810nm - 300mW pointed (module)
 650nm - 40mW pointed (module)
 650nm - 200mW cluster (diode)
 Treatment time - 0-60mins
 Laser classification - Class 3B according to IEC-60825-1

TECHNICAL DATA

Mains voltage - 100-240V / 50 or 60Hz
 Dimensions - LxBxH - 28x22x9cms
 Weight - 1.4Kgs