

ELECTRONIC TRACTION



INDOTRAC-2 is an advanced technology solid state circuitry Static / Intermittent traction unit. Moulded Aluminium gear box and metal gears gives long life to the equipment. It has silence smooth and powerful motor to pull the appropriate weight. The treatment time is controlled by digital timer and has indicators for Rest / Hold time. Separate tension control for cervical and pelvic.

FEATURES

- Imported sturdy & noise less motor
- Moulded aluminium gear box
- Patient safety remote
- Easy to operate

TECHNICAL DATA

- Mains Voltage - 220V/ 50 Hz
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H - 28 x 31 x 28 cms
- Weight - 11 Kgs

SPECIFICATIONS

- Treatment Types - Continuous and Intermittent
- Traction force - 3 - 60 Kgs (120Kgs with Doubler)
- Hold Time - 20, 40, 60, 80 sec
- Rest time - 1, 5, 10, 20 sec
- Timer - 0- 99 min. At the end of treatment time Acoustic signal and automatic reduction of traction force

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.

Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773

Web : www.hms.co.in email : sales@hms.co.in