

PREMIUM COMBO ELECTROTHERAPY



The **UNISTIM-5s** is a smart and premium combo electrotherapy. It has more than six different types of electrotherapy applications. In the development of the unit simplicity of operation and functionality were the principal consideration. All the available types of current and therapy are clearly displayed on a large LCD screen. The corresponding parameters are easily set with five soft touch keys. Pre-programmed memory with treatment applications offers protocols for more than 77 frequently occurring indications. And users can store their settings in user memory.

FEATURES

- Both medium and low frequency currents
- Russian stimulation
- 77 Pre-programs condition wise with protocols
- Cable and electrode checking mode
- User settable memory
- Soft touch Rubber keys
- Light weight and compact design
- Easy to operate
- Plastic moulded cabinet
- Large graphic LCD display
- International standards and specifications
- Advanced patient safety system
- Output mode CC / CV
- Independent intensity controls for each channel

SPECIFICATIONS

INTERFERENTIAL

- | | | |
|-------------------|---|--|
| Carrier frequency | - | 2000 Hz / 4000 Hz |
| Base frequency | - | 1 - 150 Hz (1 Hz step) |
| Sweep frequency | - | 0 - 150 Hz (1 Hz step) |
| Sweep modulation | | |
| Programs | - | 1/1S 1/5, 1/5S 6/6S |
| Therapy Modes | - | Four pole Interferential
Four pole Vector 90°
Four pole Vector 45°
Two pole pre-modulated |

RUSSIAN CURRENT

- | | | |
|----------------|---|--------------------------|
| Frequency | - | 2500 Hz |
| Ramp Rise/Fall | - | 1 - 10 Sec (1 Sec step) |
| Burst | - | 50 - 250 Hz (10 Hz step) |
| Hold ON | - | 1 - 99 Sec (1 Sec step) |
| Hold OFF | - | 1 - 99 Sec (1 Sec step) |

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in

