

COMPUTERISED VACUUM THERAPY



ACCESSORIES



DIGIVAC is an advanced Computerised Vacuum therapy unit with continuous and pulsed output. Can be integrated with UNISTIM, DIGIMED, SONOMED and INDOMED series.

FEATURES

- Continuous & Pulsed
- Low Noise & vibration
- Flexible Output cable & Electrodes
- Light Weight, User friendly
- Connect with any Model IFT
- Bio-Compatible Electrodes
- Universal Input (90-270V AC)

SPECIFICATIONS

- Vacuum Pressure - 100 to 600 m bar (Continuously Adjustable)
- Output Mode - Pulsed / Continuous
- Number of Electrodes - Four

TECHNICAL DATA

- Mains Voltage - 100–240 V / 50/60 Hz
- Dimensions - L x B x H – 30 x 20.5 x 10 cms
- Weight - 2.2 Kgs