


DIGITAL TRACTION



INDOTRAC D is advanced Cervical/Lumbar traction unit from HMS. During the development, we take care of Patient safety and Durability of the device. INDOTRAC D is fully microcontroller controller device, it has large LCD display for parameters and easy to operate.

FEATURES

- Static and Intermittent mode.
- Cervical and Lumbar traction.
- Two speed selections.
- High force and Low force can be set.
- 1 Kg step force setting.
- Wide input voltage range.
- User defined Hold, Rest and Treatment time.
- LCD display for all parameters.
- Actual force displayed in real time.
- Reliable, durable heavy-duty DC motor.
- Patient safety switch.
- Silent operation.

SPECIFICATIONS

Traction force	- 3 Kg to 60 Kg (1 Kg step)
Base force	- 0 Kg to 55 Kg (1 Kg step)
Traction mode	- Cervical and Lumbar
Treatment time	- 01 to 99 min (1 min step)
Hold time	- 01 to 99 sec (1 sec step)
Rest time	- 01 to 99 sec (1 sec step)
Traction speed	- Speed 1 (1 kg/sec) Speed 2 (3 kg/sec)

TECHNICAL SPECIFICATIONS

Input Voltage	- 100 to 240 V AC 50/60 Hz
Electrical class	- Class 1 Type B (IEC 60601-1)
Power consumption	- 20 Watts (max)

ACCESSORIES



DEVICE DIMENSIONS

Length	- 275 mm
Width	- 300 mm
Height	- 275 mm
Weight	- 9 Kgs

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.

Ph : +91 44 2476 0756, 2476 6576.

Web : www.hms.co.in email : sales@hms.co.in