

MICRO STIMULATOR WITH TENS - POCKET MODEL



DIGISTIM EMS+T is a specially designed Pocket model Stimulator with TENS unit for physiotherapy students and home visit therapists. It has four modes IG, F, SF & TENS. IG duration / SF intervals can be varied by single control. Micro-controller based circuitry helps us to make it compact and sleek. It can be operated by battery/mains.

FEATURES

Battery / Mains operated
Microcontroller based circuitry
Handy model

SPECIFICATIONS

Treatment Models - IG, F, SF
IG - 30-150 mSec.
FARADIC
Freq./P.width - 50 Hz / 0.7 mSec
Surge Duration - 1-6 sec
TENS Frequency - 2 - 125 Hz
Pulse width - 100 μ sec.
Power Source - 9 V Battery / Adaptor



First time
in
INDIA

HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.
Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773
Web : www.hms.co.in email : sales@hms.co.in