

PROGRAMABLE 4 CHANNEL TENS



DIGITENS-4s is an advanced micro controller based Four Channel TENS from HMS R & D. Unique design with manual and 20 pre-programs. Used for both fitness and therapeutic purpose. Soft touch keys for selection of programs. Built in battery with 5 hrs backup indigenous design with international standards.

FEATURES


- 4 Channel Outputs
- Micro controller based circuitry
- Manual Cum preprograms facility
- Obesity Program
- Universal Input
- Built In Re-chargeable battery with 5 hrs Back-up
- ABS Cabinet-Portable & Light Weight
- 128 / 64 Graphic LCD Display for parameter setting
- Individual intensity output display for each channel

TECHNICAL DATA

Mains Input Voltage	- 90V – 270VAC
Battery	- Maintenance free
rechargeable (Optional)	
Safety Class	- Class I, Type BF (IEC 60601 -1)
Dimensions	- L x B x H – 23 x 23 x 9.5 cms
Weight	- 1.75 Kgs

SPECIFICATIONS

Tens Frequency	- 2 – 250Hz (1 Hz Step)
Pulse width	- 20 – 250 micro seconds
No.of Channels	- Four
Output mode	- Continuous/Pulsed
Frequency Sweep	- Available
Treatment Time	- 0 – 60 min
Pre-program	- 1-15 Pain Relieving Program 16-20 Fitness Program
Display	- 128 / 64 Graphic LCD Display
Power Source	- Dual Power (Mains/Battery)



HMS Medical Systems

#9, 1st Cross Street, New Colony, Porur, Chennai - 600116.

Ph : +91 44 2476 0756, 2476 6576 Telefax : +91 44 2476 3773

Web : www.hms.co.in email : sales@hms.co.in