

4 CHANNEL TENS



INDOTENS 4ch is an advanced micro controller based Four channel TENS. Unique design with manual and 9 pre-programs for therapeutic purpose. Soft touch keys for selection of programs, Indigenous design with integrated safety features.

FEATURES

- 4 Independent channels
- Micro-controller based circuitry
- user friendly
- Universal input voltage
- Manual Cum preprograms facility
- ABS Cabinet-Portable & Light Weight.

TECHNICAL DATA

- Mains Input Voltage - 90V – 270VAC
- Safety Class - Class I, Type BF (IEC 60601 -1)
- Dimensions - L x B x H – 24 x 28 x 9.8 cms
- Weight - 1.25 Kgs

SPECIFICATIONS

- Frequency - 1 – 250Hz (1 Hz Steps)
- Pulse width - 20 – 400 micro sec
- No.of Channels - Four
- Output mode - Continuous/Pulsed
- Treatment Time - 0 – 99 min
- Pre-program - P-0 Manual Settings
P-1 to P9 Pre-Programs
- Display - 5 Digit 7Segments Display
- Frequency Sweep- Available

