



IFT +RUSSIAN + TENS + EMS COMBO

COMBOSTIM is a portable combination therapy device featuring IFT, TENS, Russian, and Muscle Stimulation with user-friendly pre-set programs.

Powered by a rechargeable battery, it delivers reliable clinical-grade performance for home, clinic, or on-the-go treatment.

FEATURES

- 4 Therapies in 1 device
- Dual independent channel output
- 40 preset programs and Manual mode
- 2.8" Touch Display
- Rechargeable Battery
- For Home and Clinical use
- Light Weight and portable (130 g)
- ABS Body, Shockproof and Heat proof



Device

SPECIFICATIONS

INTERFERENTIAL Waveform

Waveform	- Bi-phasic square Wave pulse
Carrier frequency	- 4000Hz
Base frequency	- 1 - 150 Hz
Sweep frequency	- 0 - 150 Hz
Treatment time	- 60 min
Output mode	- 4 Pole & 2 Pole

RUSSIAN

Waveform	- Bi-phasic square Wave pulse
Carrier frequency	- 2500 Hz
Output mode	- Constant/Surge
Burst Frequency	- 30 - 100Hz (10Hz step)
Ramp time	- 1 - 6 sec (1sec step)
Surge ON/OFF	- 1 - 30 sec(1sec step)
Treatment time	- 60 min

TENS

Waveform	- Bi-phasic square Wave pulse
Modes	- Continuous, Burst, PWM & PFM
Frequency	- 1 - 250 Hz(1Hz step)
Pulse width	- 50 - 400µs (10µs step)
Burst Frq	- 0.5 - 5 Hz (0.5Hz step)
Treatment time	- 60 min

EMS

Waveform	- Mono-phasic square Wave pulse
Modes	- Synchronuous, Alternate, Delay & IG pulse
Frequency	- 1 - 150 Hz(1Hz step)
Pulse width	- 50 - 700µs (10µs step)
Ramp	- 1 - 6 sec (1sec step)
ON time	- 1 - 30 sec(1sec step)
OFF time	- 0 - 60 sec (1sec step)
IG pulse	
ON time	- 10 - 100ms (10ms step)
OFF time	- 0.5, 1 & 2sec
Treatment time	- 60 min

Technical Data

Display	- 2.8" Touch display
Battery charging	- Type-C
Charging time	- 2h 30min
Usage time	- Up to 4 Hours (continuous use) (Depends on Mode and intensity)
Power source	- 3.7V/1000 mAh (Li-po)
Product Dimension	- 130 x 69 x 26 in mm
Product weight	- 130 g
Safety class	- Class II, Type BF(IEC 60601 - 1)

ACCESSORIES

